

HW-EX124

SUCTION DEVICE



Healthward

Healthier your life

HW-EX124

FEATURES AND HIGHLIGHTS

Technical parameters

 Power voltage 	AC110-220V±22V	
 Rated frequency 	50-60Hz	
Suction pump	Oil-free piston pump	
Negative pressure	≥0.08MPa(680mmHg)	
Range of negative pressure	0.02MPa-limit negative pressure value	
• Noise	≤60dB	
Suction rate	≥20L/Min	
 Storage bottle capacity 1000mL/piece, high temperature resistant glass 		
 Intermittent loading continuous operation Short operation time, maximum operation time 30min, duration 50% 		
 Classification of safety protectic Class I equipment, type B appli 		
Environmental temperature ra	nge −40~60°C	
Relative Humidity Range	≤90%	
Atmospheric pressure range	500∼1060hPa	
Weight	3.5kg	
Package size	360×270×290mm	

Technical configuration

\square Suction flexible catheter(length 2M, Φ 7×13)	1 pc
☑ Empty foot switch	1 pc
☑ Air filter	2 pcs
☑ Fuse	2 pcs
Chest suction tube	1 pc



Feature

Low noise, maintenance-free, high pumping rate, fast negative pressure rise, no oil mist pollution.

Design

Portable medical suction device, easy to move, compact and flexible. It can be used to suction out the thick sputum and thick liquid from the pharynx of obstructed patients or pediatric patients.

Structure feature

The product structure primarily consists of piston pump, vacuum gauge, negative pressure regulating valve, air filter, storage bottle, and foot pedal switch.

Working principle

Utilizing a high-efficiency piston pump, it operates without the need for oil lubrication, ensuring long-term operation and has a long service life. The structure is rational and will not generate positive pressure, absolutely guaranteeing safety.

Piston pump

Utilizing a high-efficiency piston pump, it operates without the need for oil lubrication, enabling long-term operation and has a long service life.

Shell material

Featuring an ABS injection-molded shell, it has an attractive appearance, stable and durable, and ease of cleaning.



Handle

Fixed, portable handle for easy lifting, saving time and effort.



Vacuum gauge

Using high precision vacuum pressure gauge, direct display of negative pressure value.



Negative pressure regulating valve

Any negative pressure value can be selected within the limit range for easy operation. Negative pressure adjustment range: 0.02-0.09MPa (680mmHg).



Air filter

Air filter reduce air pollution of the surrounding environment from dirty liquids.



Storage bottle

Large caliber, large-capacity liquid storage bottle, made of glass, easy to clean, with a sealing ring stopper, good sealing, easy to open. The bottle is equipped with an overflow protection device to prevent dirty liquid from flowing into the pump body.



Switch

Manual switch for low voltage control, and fully sealed, safe and reliable.


















export@health-ward.com

www.health-ward.com

