



HWA2000

ELECTRIC
OPERATION TABLE



Healthward

Healthier
your life

HWA2000

FEATURES AND HIGHLIGHTS



Technical parameters


External size (LxWxH):	2100x480x(680-930)mm
Reverse trendelenburg:	0-22°
Trendelenburg:	0-20°
Lateral tilt left/right:	0-15°
Head plate up/down:	55°/90°
Back plate up/down:	50°/15°
Leg plate up/down:	45°/90°
Leg plate outwards:	0-90°
Kidney bridge elevation:	0-100mm
Safe working load:	150kg
X-rays translucent:	NO
C-arm compatible:	NO
Multi-functional orthopedic frame:	YES

Technical configuration

<input checked="" type="checkbox"/> Bed frame	1 pc
<input checked="" type="checkbox"/> Mattress	1 set
<input checked="" type="checkbox"/> Anesthesia screen frame	1 pc
<input checked="" type="checkbox"/> Shoulder holder	2 pcs
<input checked="" type="checkbox"/> Arm holder	2 pcs
<input checked="" type="checkbox"/> Leg support	2 pcs
<input checked="" type="checkbox"/> Foot support	2 pcs
<input checked="" type="checkbox"/> Hand controller	1 pc
<input checked="" type="checkbox"/> Sliding block	10 pcs
<input checked="" type="checkbox"/> Spanner	4 pcs


Material

The base material for the table is high thickness 304 stainless steel, seamless welding by robot, with modern appearance and more stable.



Mattress

Removable memory mattress, good rebound resilience, anti-static electricity and anti-microbial.



Brake

Pedal brake, easy to move and lock.



Safety

The table is provided with an emergency stop switch on the base column, help to stop the operating table running in case of emergency situation, ensure the safety of operation.



Leg plate

Leg plate with imported gas spring for easy control.



Kidney bridge

Manual kidney bridge, more convenience.



Battery could support

Battery could support about 50-80 complete operation cycles. Working time could be last for one week according to useage and help to keep the operating table running in the event of power failure.



Armrest



Anesthesia Frame




Leg Crutches



Shoulder Support



Body Supporter



Professional functions

Lowest position of table top

680mm



Highest position of table top

930mm



Back plate up

0°-50°



Back plate down

0°-15°



Reversed trendelenburg

0°-20°



Trendelenburg

0°-22°



Right lateral tilt

0°-15°



Left lateral tilt

0°-15°



Head plate upward

0°-55°



Head plate downward

0°-90°



Leg section evagination

0°-90°



Leg plate up/down

-45°-90°



Kidney bridge elevation

0-100mm



As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.

Healthward



 17800 CASTLETON ST STE 180.CITY OF INDUSTRY CA 91748

 +1 (603) 416-3063 ·  24/7

 export@health-ward.com

 www.health-ward.com



General catalogue