

# HWA98-5

ELECTRIC OBSTETRIC BED



**Healthward** 

Healthier your life

# HWA98-5

## FEATURES AND HIGHLIGHTS













## Technical parameters

<ul><li>External size (LxW):</li></ul>	2140x1080x(750-900)mm
Height(without mattress):	580-790mm(±10mm)
Back-rest Adjustment:	0-70°
Reverse Trendelenburg:	0-6
Trendelenburg:	0-6.5
Overall height adjustment:	750mm-900mm ±10mm
Leg-rest Adjustment:	0-60°
<ul><li>Power supply:</li></ul>	AC220\
• Input power:	250VA

## Technical configuration

<b>√</b>	] Mattress	1 pc
<b>√</b>	] Basin	1 pc
<b>√</b>	] Hand controller	1 pc
_	] Handle	2pcs
_	Leg Support	2pcs
_	] Castors	4pcs
<b>√</b>	Central Locking Pedal	2pcs

# **QUALITY CONTROL**

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating. efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



#### Design

The overall design is humanized and modern, and it is an obstetric and disease comprehensive bed integrating obstetric bed, hospital bed, examination bed and operating bed.



#### **Bed material**

The bed frame is welded by 30\*60mm, ≥1.5mm high-quality rectangular steel, beautiful and burr-free, with strong bearing capacity. It can bear up to 200kg.



#### **Paintwork**

The surface of the bed body is treated with electrostatic spraying, full anti-rust treatment, anti-aging, not easy to rust, antibacterial and environmentally friendly.



#### Welding | Craftsmanship

The overall welding is adopted by the robot, and the appearance is smooth, stable and reliable.



#### **Back-rest**

The back section can be raised and lowered to assist pregnant women to sit and stand, and relieve the pressure of pregnant women during labor.



#### Leg-rest

The leg section can be lifted and lowered to assist pregnant women to give birth smoothly, and the posture can be adjusted arbitrarily, such as sitting, kneeling, lying down, and squatting.



#### Trendelenburg&Reverse Trendelenburg

Electric control is used for trendelenburg and reverse trendelenburg.



#### **Height Adjustment**

The overall height of the bed body can be adjusted to facilitate the adjustment of the height of the bed, which is more suitable for nursing staff to take care of pregnant women.





# Head & foot board European-style head and foot board, resistant to fading and aging, high strength, easy to clean, firm and strong; hidden locking switch, stable and reliable, easy to disassemble, more suitable for the concept of a family-oriented and warm delivery room.















## **Healthward**

As the circumstance of printing, the photos might have minor difference on color, tactile from real products. As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.















export@health-ward.com

www.health-ward.com

