

HVVQ5k DIALYSIS BED



Healthward

Healthier your life

HWQ5k

FEATURES AND HIGHLIGHTS







Technical parameters

| External size (LxWxH): | 1930x880x600 mm |
|--|-----------------|
| Mattress platform: | 1710x780 mm |
| Safe working load: | 200 kg |
| Back-rest adjustment: | 0-75° (±5°) |
| Knee-rest adjustment: | 0-45° (±5°) |

Technical configuration

| ☑ Linear motor | 2pcs |
|---|------|
| ☑ Aluminium alloy foldable side rail | 1set |
| ✓ PP bed ends with safe lock | 1set |
| ✓ PP bed platform, removeable, washable | 5pcs |
| ☑ 5" double side castors | 4pcs |
| ☑ IV pole | 1set |
| ✓ IV pole prevision | 4pcs |

QUALITY CONTROL

Healthward, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. Healthward adopt robot welding, automatic epoxy coating. efficient assembly line, dedicated operator, rigorous quality control, only to guarantee Healthward's commitment on quality



Design

This Bed uses unique design to meet hospital strict gremfree request. Head & Foot board and platform are all seamless, which supports effective sterilization.



Steels | welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



Plastic | Painting

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa

11 Process epoxy painting, ASTM testing anti-baterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Electric Motor

Linear motor, with UL /ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.

Accessories



Platform

Five sections PP Platform with Integrated mattress Retainer, removable for sufficient sterilization.



Castors | Brake padel

 $5^{\rm w}$ Double covered castors with TPR tire, no worn-out after running 30KM, united forming without bolts

Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500time. The double-folded large-area rubber design, horizontal and vertical two-sided rubber, allows the nurse to touch the brakes when the soles of the feet are soft and large-area plastic surface, saving effort and protecting the upper.













Safe working load

 200_{kg}



Back-rest adjustment

75°



Knee-rest adjustment

45°

Healthward

As the circumstance of printing, the photos might have minor difference on color, tactile from real products. As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.















export@health-ward.com

www.health-ward.com

