



HWA501

ELECTRIC-HYDRAULIC
OPERATING TABLE

Healthier
your life

HWA501

FEATURES AND HIGHLIGHTS



Technical parameters

External size (LxWxH):	2125x550x(600-950)mm
Reverse trendelenburg:	0-20°
Trendelenburg:	0-15°
Lateral tilt left/right:	0-15°
Head plate up/down:	55°/90°
Back plate up/down:	89°/35°
Leg plate up/down:	15°/90°
Leg plate outwards:	90°
Kidney bridge elevation:	0-120mm
Horizontal sliding:	0-320mm
Horizontal rotation:	180°

Technical configuration

<input checked="" type="checkbox"/> Bed frame	1pc
<input checked="" type="checkbox"/> Mattress	1set
<input checked="" type="checkbox"/> Anesthesia screen frame	1pc
<input checked="" type="checkbox"/> Shoulder holder with cushion	2pcs
<input checked="" type="checkbox"/> Arm holder with cushion	2pcs
<input checked="" type="checkbox"/> Leg support with cushion	2pcs
<input checked="" type="checkbox"/> Foot support with cushion	2pcs
<input checked="" type="checkbox"/> Hand controller	1pc
<input checked="" type="checkbox"/> Kidney bridge handle	1pc
<input checked="" type="checkbox"/> Battery back-up	1pc


Material

The base column and side holder are all made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.



Dual control system

All functions controlled by both hand controller or control panel on table base column.




Special functions

With table top 180° rotation function.




Battery back-up

Equipped with a built-in emergency power supply, the emergency power supply can be automatically charged, and the emergency power supply can automatically jump to the working state when the network power supply is cut off.




Brake

Electric Control system for the brake, it is easily for moving operating table.




Kidney bridge

Built-In Kidney Bridge, More Convenient And Nice Appearance.



Mattress

PU cover, memory mattress, easy cleaning, thickness of mattress: 5cm.



Anesthesia screen frame




**Leg support with cushion
Foot support with cushion**



Arm holder with cushion



Shoulder holder with cushion



Body holder with cushion



Professional functions

Lowest position of table top
600mm



Highest position of table top
950mm



Back plate up/down
89°/35°



Kidney bridge elevation
0-120mm



Reversed trendelenburg
0°-20°



Trendelenburg
0°-15°



Left lateral tilt
0°-15°



Right lateral tilt
0°-15°



Head plate upward
0°-55°



Head plate downward
0°-90°



Leg section evagination
0°-90°



Leg plate up/down
15°/90°



Horizontal sliding
0-320mm





As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.

Healthward



 17800 CASTLETON ST STE 180.CITY OF INDUSTRY CA 91748

 +1 (603) 416-3063 ·  24/7

 export@health-ward.com

 www.health-ward.com



General catalogue