## 150 D ( E PA 🎡 🦄



# Healthward `





FEATURES AND HIGHLIGHTS



### Technical parameters

External (LxW):	1800x600mm
Table height:	700-950mm
<ul> <li>Back plate up folding:</li> </ul>	0-60°
Reverse Trendelenburg:	0-25°
<ul> <li>Trendelenburg:</li> </ul>	0-16°
Load bearing:	≤150kg



### Technical configuration

🗹 Arm support	2 pcs
☑ Handle	2 pcs
☑ Leg support	2 pcs
☑ Stirrup	2 pcs
☑ Slider	4 pcs
☑ Tightener	2 pcs
☑ Hand controller	1 pc
☑ Foot switch	1 pc
☑ Stainless steel basin	1 pc







#### Material

The chassis cover is made of high quality 304 stainless steel, which is welded by the robot seamlessly. The shape is beautiful and stable.



Mattress

**Braking system** 

of the casters.

Battery

Polyurethane foam mattress, waterproof and seamless, easy cleaning and disinfection.

The pedal controls the brake function

The delivery bed is equipped with a standard battery, which could support the electric operation for 50-80

required by the operator under special circumstances such as power failure.

The delivery bed is equipped with a

stainless steel basin, which can be

quality three-stage slide rails.

stored under the bed through a high-

times, and could keep maintaining the operation of various functions



#### Hand controller

All the functions are controlled by hand controller.



#### **Control System**

Electric control the height lifting, trendelenburg and reverse trendelenburg, back plate folding up, easy and convenient operation.



#### **Pull handles**

Both sides of the delivery bed are equipped with collapsible handle, which is convenient for the puerpera to get on and off the bed, with antislip lines, to prevent slipping when the force is applied.



#### **Auxiliary Table**

The auxiliary table can be retractable and stored under the bed body, easy to use.



### Leg supports

Sewage basin

The leg support is integrated with stirrup function, and the delivery bed is equipped with a pair of stirrups.





#### **Arm Support**

The arm support could be adjusted into multiple angles, it can be used to comfortably fix the puerpera's arm at various angles.





## Healthward `

As the circumstance of printing, the photos might have minor difference on color, tactile from real products. As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.



- 17800 CASTLETON ST STE 180.CITY OF INDUSTRY CA 91748
- 🕕 +1 (603) 416-3063 · 🕓 24/7
- export@health-ward.com
- () www.health-ward.com

