

HWR058-ET

EMERGENCY TROLLEY



Healthward[®]

Healthier
your life

HWR058-ET

FEATURES AND HIGHLIGHTS

Material

The ABS material is an innovative engineering plastic. Antibacterial, moisture-proof, anti-corrosion, wear-resistant, non-toxic and tasteless, never rust. Washable, easy cleaning.



Design

The overall design and production are exquisite, the styling and color matching are reasonable, installation is simple, and the structure is stable and durable.



Drawer

Boosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



Drawer compartment

There are transparent acrylic medicine grids inside the drawer, which can be used to classify and place medicines, and the medicine grids can be freely combined.

Technical parameters

Overall Dimensions	790W X 610D X 1000H mm
Defibrillator Stand	380 X 340mm
Work Board	340 X 250mm
Cpr Board	475 X 595 X 8.5mm
Packing Size	820 X 650 X 1150 mm
Product Weight	65 Kg
Caster	5" Two Casters With Brake, Two Without Brake
Work Board Size	340 X 250mm
Bracket Size	380 X 340mm


Technical configuration

- 20L Trash Can With Lid
- Oxygen Bottle Rack
- Adjustable Iv Stand
- Medicine Grid Horizontal And Vertical
- Slide Rail

Accessories


IV stand

Adjustable aluminum alloy IV stand.



Caster

At the bottom, it uses two 5-inch silent brake, anti-winding, high-end, silent casters, allowing for flexible movement and a high load-bearing capacity.




Drawer slide

The drawer uses high-strength three-section silent damping slides. When the drawer is closed to the final distance, hydraulic pressure is used to slow down the speed and reduce the impact force.




Working table

The main tabletop can have a replaceable soft pad that is dust-proof and easy to clean. The vehicle is equipped with a hidden sliding writing board that can be used for writing.



Sealed handle

The trolley is equipped with a sealed handle, which can lock five drawers at the same time.



Oxygen bottle rack

ABS oxygen bottle bracket on the left side.



Trash can


Clamshell ABS trash can on the right side, Specification: 20L



CPR board

The ABS CPR board is installed on the back of the trolley body, which can be used during defibrillation rescue.

CPR board size: 475 * 595 * 8.5mm



As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.

Healthward



 17800 CASTLETON ST STE 180.CITY OF INDUSTRY CA 91748

 +1 (603) 416-3063 ·  24/7

 export@health-ward.com

 www.health-ward.com



General catalogue