



HWA98-3

ELECTRIC OBSTETRIC BED

Registration code: SKL-A



Healthward[™]

Healthier
your life

A98-3

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External size (LxWxH)	2120 x 1120 x (595-905) ±10mm
▪ PP back section	775x875mm
▪ PP seat section	430x880mm
▪ PP leg section	650x794mm
▪ Back-rest adjustment	0-70°
▪ Trendelenburg	0-10°

Technical configuration

<input checked="" type="checkbox"/> Leg Holder (PU material)	1pair
<input checked="" type="checkbox"/> Grab Handle	1pair
<input checked="" type="checkbox"/> Basin	1pc
<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> IV Pole Prevision	2pcs
<input checked="" type="checkbox"/> Controller	1pc
<input checked="" type="checkbox"/> Satinless Steel Guardrail	1pair
<input checked="" type="checkbox"/> 6" Castor	4pcs
<input checked="" type="checkbox"/> 80mm Mattress	1set

Material

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.

Bed board

Bed surface has PP cover, made by virgin material. The PP cover is detachable, and easy to clean.

Head&Foot Board

The bed platform and head&foot board made by virgin PP material, the head&foot board is movable.

Hand controller

Equipped with one hand controller.

Guardrail

Handrail is rotatable, when you don't use it, you can put down the guardrail.

Handle

The hand bar can be dismantled and movable.

Lithotomy rod

The leg support can be dismantled and movable.

Foot support

The foot support can be dismantled and movable.

Castor

6 inch central controlled castor for moving obstetric table.

Accessory

Equipped with stainless steel sewage basin and other various accessories.

Mattress

The mattress adopts high density memory sponge and is once formed. Mattress is movable.

Foot section

foot section can be retracted into the bed, convenient for gynecological and delivery surgical.



External Size (LxW)

2120x1120 mm



Back-rest adjustment

70°



Height adjustment

595-905 mm



Trendelenburg


10°



Foot section into the bed



CONSTANT INNOVATION




Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

Healthward

Chengdu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road ·
Wuhou District (Chengdu) · China

☎ +86 18021231901 · 🌐 24/7

✉ export@health-ward.com

🌐 www.health-ward.com



Company video



General catalogue