

HWK2k TWO FUNCTIONS MANUAL BED




Healthward`



HWK2k

FEATURES AND HIGHLIGHTS





Technical parameters

External size (LxW):	2160x980 mm
Mattress platform:	1930x830 mm
Safe working load:	200 kg
Height adjustment:	505mm
Back-rest adjustment:	0-80° (±10°)
Knee-rest adjustment:	0-35° (±10°)

Technical configuration

🗹 Manual Cranks	2pcs
☑ 5" Covered Castors	4pcs
☑ Foldable Side Rail	1set
ABS Bed Ends with Safe Lock	lset
☑ IV Pole	1set
☑ IV Pole Prevision	4pcs
🗹 Drainage Hook	4pcs
Mattress Holder	2pcs

Optional color (MOQ100pcs- color changing)



Red oak

QUALITY CONTROL



Design

The design is both modern and practical.

practical. The head and foot boards feature a beautiful and elegant shape. The use of a hollow protective wrap angle, combined with a barreltype bumper, makes the product instantly recognizable.



Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



Plastic | Painting

Environmental plastic can be sterilized at a temperature of 100°C and has a tensile strength of up to 30MPa.

11 Process epoxy painting with a thickness of 0.12mm and a brightness of 60°, which is capable of resisting impacts of up to 50kg.

Accessories



Head & foot board

The head and foot boards feature "integrated protective corners" with hollow design, which are combined with barrel-type bumper wheels that have a diameter of 60mm.

This design not only minimizes the overall space occupied by the bed but also provides protection.

Additionally, the head and foot boards are effortless to disassemble and clean.



Castors

5" covered casters.

PR material, walks 30KM with no wear and tear. Precision-designed antiwinding shell with no screws.

Pass Dynamic test: withstand a weight of 120kg, walk 30KM, and pass through 500 obstacles without any deformation or damage.

Platform

Hand Cranks

100,000 uses.

Two retractable hand cranks to control backrest and kneerest. The hand cranks remain intact after

Die punching bed board features a robot-welded surface that has been polished to perfection, leaving no welding scars. Bearing 200kg.





Mattress Retainer

The dual-function mattress retainer serves to not only prevents scratches on the bed frame's paint but also retrains the mattress to its place, and prevents mattress sliding.



IV Pole

The infusion pole is made of stainless steel and has four plastic hooks that can bear up to 15 kg of weight.





Safe working load





Back-rest adjustment





Knee-rest adjustment



Healthward`

As the circumstance of printing, the photos might have minor difference on color, tactile from real products. As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.



♥ 17800 CASTLETON ST STE 180.CITY OF INDUSTRY CA 91748

🕕 +1 (603) 416-3063 · 🕓 24/7

export@health-ward.com

www.health-ward.com

