ISO :: D ( E FDA 🎬 🦄

# K3k THREE FUNCTIONS MANUAL BED

Registration code: SK-CD1



## Healthward`



## K3k FEATURES AND HIGHLIGHTS





#### Technical parameters

External size (LxW):	2105x980 mm
Mattress platform:	1925x820 mm
<ul> <li>Safe working load:</li> </ul>	200 kg
<ul> <li>Height adjustment:</li> </ul>	400-710mm
<ul> <li>Back-rest adjustment:</li> </ul>	0-75° (±10°)
Knee-rest adjustment:	0-40° (±10°)

### Technical configuration

3 pcs
4 pcs
1 set
1set
1set
4pcs
4pcs
2pcs

#### Optional color (MOQ100pcs- color changing)



Red oak

HEALTHWARD

# **QUALITY CONTROL**



#### Design

The design is both modern and practical. The head and foot boards feature a beautiful and elegant shape. The use of a hollow protective wrap angle, combined with a barreltype bumper, makes the product instantly recognizable.



**Steels | welding** Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



#### Plastic | Painting

Environmental plastic can be sterilized at a temperature of 100°C and has a tensile strength of up to 30MPa.

11 Process epoxy painting with a thickness of 0.12mm and a brightness of 60°, which is capable of resisting impacts of up to 50kg.

### **Accessories**



#### Head & foot board

The head and foot boards feature "integrated protective corners" with hollow design, which are combined with barrel-type bumper wheels that have a diameter of 60mm.

This design not only minimizes the overall space occupied by the bed but also provides protection.

Additionally, the head and foot boards are effortless to disassemble and clean.



#### Castors

5" covered casters.

PR material, walks 30KM with no wear and tear. Precision-designed antiwinding shell with no screws.

Pass Dynamic test: withstand a weight of 120kg, walk 30KM, and pass through 500 obstacles without any deformation or damage.

#### Platform

Die punching bed board features a robot-welded surface that has been polished to perfection, leaving no welding scars. Bearing 200kg.





Three retractable hand cranks to control back-rest, knee-rest and height. The hand cranks remain intact after 100,000 uses.





Siderail The guardrail is made of six-pole folding aluminum alloy, with stainless steel poles and a plastic handle.

#### **Mattress Retainer**

The dual-function mattress retainer serves to not only prevents scratches on the bed frame's paint but also retrains the mattress to its place, and prevents mattress sliding.



The infusion pole is made of stainless steel and has four plastic hooks that can bear up to 15 kg of weight.





Back-rest adjustment





Knee-rest adjustment





Height adjustment









Design







Manufacturing



Maintenance





#### USA

Saikang Medical Technology, Inc. 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

#### United Kingdom

UK Healthward International Ltd. Dept. 111, 196 High Road Wood Green, N22 8HH, London. Tel.: +44 2032399738

#### Russia

Saikang Medical Russia LLC Truda 174, 454080 Chelyabinsk. Tel.: +7 (351) 214 5566 / 8811 / 1516

#### East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel.: +254 733 704 690

#### Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497

## Healthward`

Chengdu Saikang Medical Equipment Co., Ltd. STOCK CODE: 870098

### 📀 No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road ·

- □ +86 18021231901 · S 24/7
- 💌 export@health-ward.com
- () www.health-ward.com







General catalogue