# HWV2k MANUAL BED




### Healthward`

ISO :: C ( FDA 🎬 🦄



# HWV2k

#### FEATURES AND HIGHLIGHTS





#### Technical parameters

<ul> <li>External size (LxW):</li> </ul>	2130x995 mm
Mattress platform:	1925x900 mm
Safe working load:	200 kg
Height adjustment:	500mm
<ul> <li>Back-rest adjustment:</li> </ul>	0-70° (±10°)
Knee-rest adjustment:	0-38°(±10°)

#### Technical configuration

🗹 Manual crank	2 pcs
☑ 5'' covered castors	4 pcs
☑ Foldable side rail	1 set
☑ PP bed ends with safe lock with Φ65 bumper wheels	1 set
PP bed platform with integrated mattress retainer	1 set
☑ IV pole 4 hooks	1 set
☑ IV pole prevision	4 pcs
🗹 Drainage hook	2 pcs

#### Optional

Auto-regression function (back section)
Utility shelf
Trapeze bracket/pole

#### Optional color (MOQ100pcs- color changing)



Red oak

## **QUALITY CONTROL**

Healthward, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. Healthward adopt robot welding, automatic epoxy coating. efficient assembly line, dedicated operator, rigorous quality control, only to guarantee Healthward's commitment on quality



Design

The unique germ-free design meets hospitals' strict infection control requirment. Seamless PP head and foot board and platform ensure the instant and sufficient sterilization



**Steels | welding** Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



Plastic | Painting Environmental plastic can be sterilized in 100°C, tensile upto 30MPa. 11 Process epoxy painting, ASTM testing anti-baterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.

#### Accessories



#### Platform

Integrated PP bed platform, adopting advanced manufacturing technology. 4-5 section platform, enables independent installation and easy maintenance.

 $\ensuremath{\mathsf{Enlarge}}$  the seat area when back part is raised, providing more comforts for the patients.

Reduce the depression of the sponge at the seat panel.



#### Castors

5" covered casters. TPR tire no worn out after running 10 km. Anti-winding hard shell cover,united forming without bolts. Pass dynamic test: pass obstacles 40000 times, under the standard 200kg load of the whole bed.



LG Chem fresh PP material, Plug in type. Equipped with the bumpers wheel, diameter 65mm, preventing collision.





Platform integrated mattress retainer at the both sides of the back and leg section, with a boarder of 33mm height.









IV Pole SS IV pole, bearing 15kg.





Safe working load





Back-rest adjustment



Knee-rest adjustment



### Healthward`

As the circumstance of printing, the photos might have minor difference on color, tactile from real products. As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.



♥ 17800 CASTLETON ST STE 180.CITY OF INDUSTRY CA 91748

🕕 +1 (603) 416-3063 · 🕓 24/7

export@health-ward.com

🌐 www.health-ward.com

