ISO 🔅 C E FDA 🎆 🦄

HWH048 OXYGEN BOTTLE TROLLEY



Healthward`



HWH048

FEATURES AND HIGHLIGHTS

Technical parameters

Size (LxWxH)

1435×360×1040 mm

Technical configuration

☑ 6" TPP Rubber Rastor	2pcs
☑ 3" Covered Castors	2pcs

Design

U-shaped groove design, the radian meets the requirements of the oxygen bottle for safe transportation.





Size: $1435 \times 360 \times 1040$ mm, Can load 40L oxygen cylinders, the body is light and easy to move.

Handle

Pushing handle, easy to move the trolley.

6" TPP Rubber Rastor

6" TPP rubber castor, Strong, durable, labor-saving, safty in use.







Steels | Welding Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth wedlding.



Galvanized Chain 4mm galvanized chain, operational safety.



Thickening Plate Cold rolled steel plate, effective compression, bearing weight, not easy to damage, long service life.

3" Covered Castors 3" Covered Castors, each with brakes, easy to install, stable and smooth.







Healthward `

As the circumstance of printing, the photos might have minor difference on color, tactile from real products. As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.



- 17800 CASTLETON ST STE 180.CITY OF INDUSTRY CA 91748
- 🕕 +1 (603) 416-3063 · 🕓 24/7
- export@health-ward.com
- () www.health-ward.com

